



## **1. Product Name**

TheraGym® Seesaw Glider & Swing

## **2. Product Code**

21550

## **3. Colour**

As shown

## **4. Brief Description**

Three swings in one! This multipurpose vinyl-covered swing can be used in a variety of ways: as a traditional glider; as a seesaw; or on an incline.

## **5. Contents**

- 1 x See Saw Swing
- 2 x Rope assemblies
- 6 x Locking Safety Snaps

## **6. Snoezelen® Stimulations**

- Vestibular

## **7. Best Use**

Use from a safe and appropriate suspension system.

## **8. Compatible ROMPA® Products**

21928 – Hanging Bracket

21559 – Safety Rotational Device

21570 – Safety Snap

See [www.rompa.com](http://www.rompa.com) for mats and a large range of Sensory Integration products

## **9. Starting Up**

- There are no short-cuts when setting up and using Sensory Integration equipment. Please keep these instructions in a safe place and refer to them before each use.

## **10. Detailed Description**

Three swings in one! This multipurpose vinyl-covered swing can be used in a variety of ways: as a traditional glider; as a seesaw; or on an incline, providing significant versatility in the therapy space. With three points of attachment for multiple functions, simply adjust the height of the rungs by knotting the rope. Comes complete with 6 safety snaps and a 182cm long rope. Working load: 180kg (400 lbs)

## **11. Safety**

- Always supervise the use of this product
- Always use the Rotational Device to help prevent wear of your swing and the suspension points.
- Do not exceed the working load: 180kg (400 lbs).
- Sensory Integration products are intended to be used by professionals.
- Use safety mats on the floor, as well as on the walls as necessary.
- Use a safety helmet, and any other appropriate safety measures if required for your client's needs.
- Use only steel safety snaps as aluminium safety snaps wear more quickly.
- For indoor use only.
- Check all hardware and stitching before each use to check these are intact.

This product is not a toy. Any activity involving motion, height, or rotation may cause serious injury. This product is intended for use under the supervision of trained professionals.

## 12. Technical Specification

Working load: 180kg (400 lbs)

Size:                    Size: 122cmL x 43cmW

## 13. Installation

Allow plenty of unobstructed space around your swing. Suspended equipment needs a minimum ceiling height of 2.5m (8 feet). Suspended equipment should be positioned a clear distance of 1.8m (6 feet) from walls and obstacles. This equipment is for indoor use only, in dry environments.



## Safety Snaps

Our Safety Snaps and the included Carabiner can be locked.

1. Turn the locking nut to open the snap (to slide this nut down to allow the snap to be opened). Turn the nut in the opposite direction to slide it up and over the opening to lock it.



2. The locking nut of the Safety Snap should not move from side to side at any great angle. If it does, discontinue use immediately.

Attach the rope assemblies to the swing using the safety snaps and the loops on the swing. There are various options for this, depending on the movement required.

Attach the top of each rope assembly to an appropriate suspension system.

Attaching the safety snaps and ropes to each corner will allow a more stable movement.



Attaching the safety snaps to the middle loops on the swing will allow a see saw movement.



The rope assemblies can also be attached to the loops at the end and in the centre.



## 14. Care and Maintenance

1. All the parts of your suspended products should be checked on a regular basis (at least monthly, and more if the swing is used often). Use a maintenance sheet (such as the example one provided) to record these checks.
2. Check all the parts for signs of wear and tear before each use and replace these as appropriate.
3. Use and store this product indoors, in a safe place. Do not let the product become wet.
4. Dust or wipe the components clean with a damp cloth. Do not immerse in water. Allow to dry before use.

## 15. Troubleshooting

*Not applicable*

Further copies of these instructions can be downloaded at [www.rompa.com](http://www.rompa.com)

We hope you find our instructions invaluable. If you have any suggestions for improving them further your comments will be greatly received – please contact us at [producthelp@rompa.com](mailto:producthelp@rompa.com)

